For Flu A/B Antigen Test
1. POSITIVNĚ:
1.1 Flu A Positive:
The presence of two lines as control line (C) and FluA test line within the result window indicates a positive result for Influenza A viral antigen.
1.2 Flu B Positive:
The presence of two lines as control line (C) and FluB test line within the result window indicates a positive result for Influenza B viral antigen.
1.3 Flu A+B Positive:
The presence of three lines as control line (C), FluA test line and FluB test line within the result window indicates a positive result for Influenza A and Influenza B viral antigen.
2. NEGATIVĚ:
The presence of only control band (C) within the result window indicates a negative result.
3. INVALID:
If the control band (C) is not visible within the result window after performing the test, the result is considered invalid. Some causes of invalid results are because of not following the directions correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be re-tested using a new test. If the problem persists, please stop using the product and contact with the test distributor.
For COVID-19/RSV Antigen Test
1. POSITIVNĚ:
1.1 RSV Positive:
The presence of two lines as control line (C) and RSV test line within the result window indicates a positive result for RSV viral antigen.
1.2 COVID-19 Positive:
The presence of two lines as control line (C) and CoV test line within the result window indicates a positive result for COVID-19 viral antigen. If tested COVID-19 positive, users should not take any decision of medical relevance without first consulting a medical practitioner.
1.3 COVID-19+RSV Positive:
The presence of three lines as control line (C), RSV test line and CoV test line within the result window indicates a positive result for RSV and COVID-19 viral antigen.
2. NEGATIVĚ:
The presence of only control band (C) within the result window indicates a negative result.
3. INVALID:
If the control band (C) is not visible within the result window after performing the test, the result is considered invalid. Some causes of invalid results are because of not following the directions correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be re-tested using a new test. If the problem persists, please stop using the product and contact with the test distributor.
Note:

- Whether the test result is positive or negative, the used test cassettes, swabs and collection tubes should be put into the collection bag and then disposed of according to local regulations.
- When the test results are positive (regardless of which pathogen tested positive), you are likely infected with Flu A/B, RSV or COVID-19. In case of positive result, you should take appropriate precautions and seek timely help from a healthcare provider. Individuals who test positive with COVID-19 virus antigen should follow local regulations and seek follow-up care with their physician or healthcare provider as additional testing may be necessary.
- When the test result is negative, it means you may not be infected with Flu A/B, RSV or COVID-19. However, there is possibility that a false negative result appears. If you are tested negative but still with symptoms present or you were ever in contact with a suspect, it is necessary to repeat the test after at least 24 hours and no more than 48 hours, or use PCR reagents, or go to the hospital to confirm whether infected with FluA/FluB, RSV or COVID-19.

LIMITATION OF THE TEST

- This test kit is only used for *in vitro* diagnosis.
- This test kit is only used for qualitative detection and cannot indicate the level of antigens in the specimen.
- This test is not a substitute for a medical consultation and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions.
- Failure to follow the instructions for sample collection and testing will lead to erroneous results, and in this case the results are invalid.
- If the antigen content in the sample is below the detection limit of the product, a false negative result will appear.
- If the test result is negative, but the symptoms still exist, you need to contact your doctor for confirmation.
- A negative test result may occur if the specimen is collected, extracted or transported improperly.
- A negative test result does not rule out the possibility of infection and will not set you free from the local rules to control COVID-19 spread (e.g. contact restrictions and protective measures).
- A positive test result cannot exclude co-infection with other pathogens.

PERFORMANCE CHARACTERISTICS

1. Limit of Detection (Analytical Sensitivity)
The LoD of Influenza A for Influenza A/B+COVID-19/RSV Combo Ag Test was 1.5x10⁴ TCID₅₀/mL, Influenza B for this kit was 1.5x10⁵ TCID₅₀/mL, COVID-19 for this kit was 200 TCID₅₀/mL and RSV for this kit was 1.0x10⁴ TCID₅₀/mL.

2. High Dose Hook Effect
No high dose hook effect was observed with up to 4.8 x 10⁷ TCID₅₀/mL of Flu A virus, or up to 4.8 x 10⁷ TCID₅₀/mL of Flu B virus, or up to 1.8 x 10⁶ TCID₅₀/mL of SARS-CoV-2, or up to 3.2 x 10⁷ TCID₅₀/mL of RSV with the Influenza A/B+COVID-19/RSV Combo Ag Test.

3. Analytical specificity
For Flu A/B Antigen Test:
3.1 Cross-reactivity
The Flu A/B Antigen Test was evaluated with a total of 36 bacterial and viral isolates. Bacterial isolates were evaluated at a concentration between 10³ and 10⁸ org/mL. Viral isolates were evaluated at a concentration between 10³ and 10⁸ TCID₅₀/mL. Adenovirus 18 and Parainfluenza virus 3 were tested at 10⁷ TCID₅₀/mL. None of the organisms or viruses listed below gave a positive result in the Flu A/B Antigen Test.

Potential Cross-Reactant		
Virus	Human Adenovirus B	Parainfluenza virus 3
	Human Rhinovirus 2	Human respiratory syncytial virus B
	Human Adenovirus C	SARS-COV-2
	Human Rhinovirus 14	Acinetobacter calcoaceticus
	Human Adenovirus type 10	Bacteroides fragilis
	Human Rhinovirus 16	Neisseria gonorrhoeae
	Adenovirus type 18	Neisseria meningitidis
	Measles	Pseudomonas aeruginosa
	Human Coronavirus OC43	Staphylococcus aureus
	Mumps	Streptococcus pneumoniae
	Human Cocksackievirus A9	Streptococcus sanguis
	Sendai virus	Proteus vulgaris
Bacteria	Coxsackievirus B5	Streptococcus sp. Gp.B
	Parainfluenza virus 2	Streptococcus sp.Gp.C
	Human herpesvirus 2	Streptococcus sp. Gp.G
	Influenza B	Mycobacterium tuberculosis
	Influenza A	Mycoplasma orale
	Human respiratory syncytial virus A	Poolled human nasal wash

3.2 Interfering Substances
Whole blood, and several over-the-counter (OTC) products and common chemicals were evaluated and did not interfere with the Flu A/B Antigen Test at the levels tested: whole blood (2.5%), three OTC mouthwashes (25%), three OTC nasal sprays (10%), 4-Acetamidophenol (10 mg/mL), Acetylsalicylic Acid (20 mg/mL), Chlorpheniramine (5 mg/mL), Dextromethorphan (10 mg/mL), Diphenhydramine (5 mg/mL), Ephedrine (20 mg/mL), Guaiaacol glyceryl ether (20 mg/mL), Oxymetazoline (10 mg/mL), Phenylephrine (100 mg/mL), and Phenylpropanolamine (20 mg/mL).

For COVID-19/RSV Antigen Test:
3.3 Cross-reactivity
Results demonstrated that the COVID-19/RSV Antigen Test has no significant cross-reactivity with the organisms or viruses listed below:

Potential Cross-Reactant		
Virus	Adenovirus	Respiratory Syncytial Virus B
	Human metapneumovirus (hMPV)	Influenza A H7N9
	Rhinovirus	Influenza B Guangdong/120/00
	Enterovirus/Coxsackievirus B4	Influenza B Yamagata
	Human coronavirus OC43	MERS-CoV
	Human coronavirus 229E	Bordetella pertussis
	Human coronavirus NL63	Chlamydia pneumoniae
	Human coronavirus HKU1	Haemophilus influenzae
	Human parainfluenza virus 1	Legionella pneumophila
	Human parainfluenza virus 2	Mycoplasma pneumoniae
	Human parainfluenza virus 3	Streptococcus pneumoniae
	Human parainfluenza virus 4	Streptococcus pyogenes (group A)
Bacteria	Influenza A H3N2	Mycobacterium tuberculosis
	Influenza A H1N1	Staphylococcus aureus
	Influenza A H5N1	Staphylococcus epidermidis
	SARS-CoV-2	Poolled human nasal wash
	Respiratory Syncytial Virus A	Yeast
		Candida albicans

3.4 Interfering Substances
Whole blood, Mucin, and several over-the-counter (OTC) products and common chemicals were evaluated and did not interfere with the COVID-19/RSV Antigen Test at the levels tested: whole blood (2.5%), Mucin (2%), Phenylephrine (15%), Sodium Chloride (5%), Cromolyn (15%), Oxymetazoline(15%), Fluconazole (5%), Benzocaine (0.15%), Veratramine (20%), Zincum gluconium (i.e., Zicam) (5%), Alkalol (10%), Fluticasone Propionate (5%), Phenol (15%), Tamiflu (Osetamivir phosphate) (0.5%), Mupirocin (0.25%), Tobramycin (0.0004%).

4. Clinical Performance
Clinical performance characteristics of the CorDx Influenza A/B+COVID-19/RSV Combo Ag Test was evaluated in the clinical studies. A total of 452 symptomatic or asymptomatic subjects were enrolled for the clinical study of the Influenza A/B and RSV test. And a study on 560 symptomatic or asymptomatic suspects of COVID-19 was conducted. The test results are as follows:
For FLU A antigen detection, the positive coincidence rate is 100.00%, the negative coincidence rate is 99.34%, the total coincidence rate is 99.43%.
For FLU B antigen detection, the positive coincidence rate is 96.00%, the negative coincidence rate is 99.67%, the total coincidence rate is 99.15%.
For RSV antigen detection, the positive coincidence rate is 98.98%, the negative coincidence rate is 99.21%, the total coincidence rate is 99.14%.
For COVID-19 antigen detection, the positive coincidence rate is 89.09%, which is 98.67% when Ct ≤ 25, 93.00% when Ct ≤ 30, and the negative coincidence rate is 100.00%, the total coincidence rate is 97.86%.

INDEX OF SYMBOLS			
	Do not re-use		Batch code
	In vitro diagnostic medical device		Use-by date
	Store at 2~30°C		Consult instructions for use
	Authorized representative in the European Union		Manufacturer
	Do not use if package is damaged and consult instructions for use		Catalogue number
	Keep away from sunlight		Keep dry
	Contains sufficient for <n>tests		CE Mark

CorDx, Inc.
9540 Waples St Unit C, San Diego, CA 92121
Manufacturing site: Core Technology Co., Ltd.
Room 100, C Building, No.29 Life Park Rd.,
Changping District, Beijing 102206, P.R. China

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